

Smart Handheld Full- HD Digital Video Colposcope



See what you have been missing before....

Portable COLpro222DX

Quality Standards:

CE

ISO 13485:2016

ISO 9001:2015

CDSCO Certified

EU-MDR Approved

Revolutionizing Cervical Screening – A compact, portable and user-friendly video colposcope designed for modern gynecological practice. Ideal for clinics, outreach programs and telemedicine use. Built by **Divine Meditech** under the **Make in India** initiative, this device bridges the gap between affordability, accessibility and advanced diagnostics.

Key Features

◆ High-Definition Imaging

Sharp, detailed visuals for accurate cervical examinations.

◆ 10X Magnification

Get closer to every detail with precision magnification.

◆ OIS & PDAF Technology

Steady, crystal-clear images even during handheld use.

◆ Portable & Ergonomic Design

Lightweight, cordless design with long battery life.

◆ Seamless Documentation

Capture, store and integrate with EMR systems effortlessly.

◆ Single Sitting Approach

Real-time examination and reporting in one visit.

◆ Acetone Response Timer

Built-in acetone response timer for consistent results.

◆ Built in LED Light Source

Bright, uniform light for enhanced visibility and accurate colors.

◆ Green Filter Option

Enhances visualization of blood vessels and tissue contrast.



Choose the proMIS Smart Colposcope for:

- **Easy Operation** – User-friendly software for data, images, videos & EMR integration
- **Faster Workflow** – Minimize paperwork, maximize efficiency
- **Affordable Excellence** – Premium performance at a lower cost than traditional colposcopes
- **Patient-Centric** – Clear visuals to improve patient understanding and engagement

Image, Video & Patient Data Management Software

The software is not only a powerful IT tool in the hands of busy doctors to manage images, capture videos, document findings and patient details but also provides them an opportunity to make use of referral management, electronic medical records and telemedicine.

Key Features:

- **Password-protected login** for patient data privacy
- **Live video streaming** + high-pixel image capture
- **Compare current and past patient images** side by side
- **Cloud-based storage** with teleconsultation support
- Ready for future application upgrades
- **Searchable patient database** for quick record retrieval
- **Annotation tools** to highlight areas of concern on images
- **Multi-format diagnosis reports (1–6 images)** – printable or shareable
- **Offline mode** for capturing and storing data without internet connectivity
- **PDF Report Saving & Printing** – save the report in PDF format and take printout
- **Easy Report Sharing** – share the report via WhatsApp or compatible email provider
- Zoom in and zoom out

Smart Connectivity & System Integration

- Secure cloud-based clinical record system
- HL7 compliant & FHIR-enabled
- API-ready – easy interface with hospital or clinic management systems
- Integrated teleconsultation support
- Designed for both basic and advanced healthcare environments



"A game-changer for cervical health — Smart Handheld Full HD Colposcope
Precision Imaging, Anytime, Anywhere."

Specifications

Portable COLpro222DX

Camera Sensor	:	High-resolution CMOS (3088 × 1440 px)
Operating System	:	Android 7
Processor	:	Octa-Core (Snapdragon 8 Gen 1 or Exynos 2200)
Display	:	Touch screen 6.8" Full HD AMOLED
Magnification	:	10X
Ratio	:	20:9
Battery	:	4700 mAh
Auto-focus	:	Phase Detection Auto Focus (PDAF)
Image Format:	:	PNG, JPEG
Connectivity	:	EMR integration, telemedicine support

Standard Accessories

Charging Data Cable	1 Pc.
Power Adaptor	1 Pc.
Ring Attachment	1 Pc.
Carry Case	1 Pc.



Contact us:-

Divine Meditech Private Limited

Manufacturing Unit-I

G-197, Sector-63, Noida Gautam Buddha
Nagar, Uttar Pradesh - 201309 (India)

☎ 9818396505

✉ info@divinemeditech.in

🌐 www.divinemeditech.in



Manufacturing Unit-II

Plot No.: 18, Phase III, Industrial Area,
Tahliwal Dist. Una, H.P. - 174301 (India)

☎ 9313033505

