



# HIGH FLOW CO2 INSUFFLATOR DVFLOW-30-M



Years of  
proven  
experience



Discover more of our product at:  
[www.divinemeditech.in](http://www.divinemeditech.in)



# CO2 INSUFFLATOR 30 LTR. METAL

## DVFLOW-30-M

The units are completely microprocessor controlled which guarantees stable pressure balance without any disruption of gas flow. The units are very simple to use and are showing the data of the pre-selected pressure (mmHg); the pre-selected Flow Lt/min and the volume of gas used by LED displays

### PRODUCT FEATURES

- completely microprocessor controlled
- Automatic flow rate management according to leakages
- Stable pressure balance
- Easy to use
- Display in Real Time
- Maximum flow – 30 ltr/min



### TECHNICAL DATA

|                       |                                              |
|-----------------------|----------------------------------------------|
| Insufflation Pressure | : 4-30 mmHg                                  |
| Flow Rate             | : 1 to 30 Lt/ min                            |
| Power Supply          | : 100 / 240 VAC, 50-60 Hz                    |
| Gas                   | : CO <sub>2</sub> , carbon dioxide (Medical) |
| Dimensions            | : 310mm x 120mm x 310mm                      |
| Weight                | : 6.5 kg                                     |
| Class                 | : Class 1 type                               |

### Accessories

- Main Power Cord
- Silicone Patient Tube with Connectors
- High Pressure Hose

### Contact Us:

**Divine Meditech Private Limited**



#### Corporate Offices:-

G-197, Sector-63, Noida, Gautam Buddha Nagar, Uttar Pradesh- 201309

☎ 91-9818396505 ✉ info@divinemeditech.in  
🌐 www.divinemeditech.in



#### Manufacturing Unit:-

Plot No.: 18, Phase III, Industrial Area, Tahliwal, Dist. Una, H.P. - 174301 (India)

☎ 91-9818396505 ✉ info@divinemeditech.in  
🌐 www.divinemeditech.in