

UROLOGY

We Support



DIVINE MEDITECH PVT. LTD.



Manufacturer:

Plot No.: 18, Phase III, Industrial Area, Tahliwal
District Una, Himachal Pradesh – 174301
India

www.divinemeditech.com

info@divinemeditech.com

Delhi Office:

F-3, Building No. 9, A-Block Market,
Surajmal Vihar, Delhi – 110092
India

+91-11-32902571

ABOUT US:

DMPL is an authentic Indian manufacturer specialized in the manufacturing of high end medical devices for various endoscopy applications. We support *“Make in India”* campaign initiated by Govt. of India & aim to be a leading Indian manufacturer in our field.

We are strongly supported by our R & D Teams in India and abroad. Quality is kept intact by adhering to International quality management systems and good manufacturing practices.

DMPL has been fairly successful in satisfying the customer needs by innovations and providing them the products they actually need to excel in their clinical applications.

For the distributors, **DMPL** has succeeded in providing marketable, high quality, cost competitive, reliable and unique specification products which are difficult to be compared by competition. Therefore the channel partners and distributors can provide excellent products to their markets at fair margins.

DMPL products are being used by leading medical colleges, universities, corporate hospitals, research institutes and medical professionals in India and abroad. Delighted customers in hospitals as well as happy distributors are a testimony to our success.

DMPL has the capacity to produce the medical devices according to the requirements of the customers and plans to concentrate on major distributors and OEM partners for further business expansion.



OEM:

“Where Outsourcing comes full circle”

As an original medical device manufacturer **DMPL** provides traditional high quality contract manufacturing services whereupon our skilled engineers can develop and design our products according to your specifications, name, logo and design from start to finish.

This facility of OEM manufacturing is especially useful to those companies who do not wish to build their own manufacturing capability or wish to outsource the manufacture of their products.

Some reputed Companies / brands are already availing the OEM capability offered by **DMPL** under OEM contract.

For OEM / ODM queries,
please contact us at:
info@divinemeditech.com

INDEX

S.NO	PRODUCT NAME	PAGE NO
01	Cysto-Urethroscope sheaths with Obturator	3
02	Cystoscope Bridge with Instrument Channels	4
03	Resectoscopes / HF Cord	5
04	Electrodes	6
05	Flexible Electrodes / Toomey Syringe	7
06	Evacuator	8
07	Urethrotome Sheath / Stone Crushing Forcep	9
08	Bladder Stone Element / Otis Urethrotome	10
09	Mauermayer Punch	11
10	Forceps for Cysto-Urethroscopy	12
11	Forcep for Nephroscopy	13
12	Accessories	14

*Color, specification, shape may vary due to ongoing product development

*Due to space constraints, pictures of some products are not included.


CYSTO-URETHROSCOPE SHEATHS WITH OBTURATOR



S.NO.	DESCRIPTION	ART NO.	SIZE
01	Cysto-Urethroscope sheath with Obturator	DM-CYS-200-01	17 Ch. / Fr.
02	Cysto- Urethroscope sheath with Obturator	DM-CYS-200-02	19 Ch. / Fr.
03	Cysto- Urethroscope sheath with Obturator	DM-CYS-200-03	20 Ch. / Fr.
04	Cysto- Urethroscope sheath with Obturator	DM-CYS-200-04	22 Ch. / Fr.
05	Cysto- Urethroscope sheath with Obturator	DM-CYS-200-05	25 Ch. / Fr.


**Telescope excluded*

Visual Obturator



S.NO.	DESCRIPTION	ART NO.	DIAMETER
01	Visual Obturator	DM-VO-200-01	17 Ch. / Fr.
02	Visual Obturator	DM-VO-200-02	19 Ch. / Fr.
03	Visual Obturator	DM-VO-200-03	20 Ch. / Fr.
04	Visual Obturator	DM-VO-200-04	22 Ch. / Fr.
05	Visual Obturator	DM-VO-200-05	25 Ch. / Fr.

CYSTOSCOPE BRIDGE WITH INSTRUMENT CHANNELS




S.NO.	DESCRIPTION	ART NO.
01	Cystoscopy Bridge Single Channel	DM-CBS-200-06
02	Cystoscopy Bridge Double Channel	DM-CBD-200-07

WORKING ELEMENT



S.NO.	DESCRIPTION	ART NO.
01	Working Element, PASSIVE	DM-WP-200-08
02	Working Element, ACTIVE	DM-WA-200-09

RESECTOSCOPES



S.NO.	DESCRIPTION	ART NO.
01	Resectoscopy Inner / Outer sheath with Obturator	DM-IO-200-10
02	Resectoscopy visual Obturator	DM-VO-200-11
03	Resectoscopy deflecting Obturator	DM-DO-200-12

HF CORD




S.NO.	DESCRIPTION	ART NO.	WORKING LENGTH
01	HF Cord	DM-201-016	Standard

ELECTRODES



S.NO.	DESCRIPTION	ART NO.
01	Loop Electrode	DM-LE-200-13
02	Ball Electrode	DM-BE-200-14
03	Collins knife Electrode	DM-CE-200-15
04	Roller Ball Electrode	DM-RB-200-16
05	Vapor Electrode	DM-VE-200-17



FLEXIBLE ELECTRODES



S.NO.	DESCRIPTION	ART NO.	DIAMETER
01	Flexible Ball Electrode	DM-FB-200-18	5 Fr.
02	Flexible Needle Electrode	DM-FN-200-19	5 Fr. 7 Fr.
03	Flexible Bugbee Electrode	DM-FE-200-20	5 Fr. 7 Fr.



TOOMEY SYRINGE



S.NO.	DESCRIPTION	ART NO.	WORKING LENGTH
01	Toomey Syringe	DM-200-100	75 ml 100 ml

EVACUATOR




S.NO.	DESCRIPTION	ART NO.
01	Bard Evacuator	DM-BE-200-22




S.NO.	DESCRIPTION	ART NO.
01	Bottle Type Evacuator	DM-BT-200-23

URETHROTOME SHEATH




S.NO.	DESCRIPTION	ART NO.
01	Urethrotome Sheath 21 Fr. with Obturator	DM-US-200-21
02	Urethrotome Bridge with one working channel 5 Fr.	DM-UB-200-22
03	Guide Tube half round	DM-GT-200-23

STONE CRUSHING FORCEP




S.NO.	DESCRIPTION	ART NO.
01	Stone Crushing Forcep	DM-CF-200-28

BLADDER STONE ELEMENT




S.NO.	DESCRIPTION	ART NO.
01	Bladder Stone Element	DM-BS-200-24
02	Visualizing Obturator	DM-VO-200-26

OTIS URETHROTOME




S.NO.	DESCRIPTION	ART NO.
01	Otis Urethrotome for adults, complete	DM-OU-200-25

MAUERMAYER PUNCH





S.NO.	DESCRIPTION	ART NO.	TYPE
01	Mauermayer Punch Complete	DM-MP-200-26	Karl Storz compatible Telescope
02	Mauermayer Punch Complete	DM-MP-200-27	Richard Wolf compatible Telescope


FORCEPS FOR CYSTO-URETHEROSCOPY





S.NO.	DESCRIPTION	ART NO.	DIAMETER	WORKING LENGTH
01	Flexible Alligator Forcep	DM-A5-300-40 DM-A7-300-40	5 Fr. 7 Fr.	40 cm 40 cm
02	Semi Rigid alligator Forcep	DM-R-A5-301-40 DM-R-A7-301-40	5 Fr. 7 Fr.	40 cm 40 cm
03	Flexible Biopsy Forcep	DM-B5-302-40 DM-B7-302-40	5 Fr. 7 Fr.	40 cm 40 cm
04	Semi Rigid Biopsy Forcep	DM-R-B5-303-40 DM-R-B7-303-40	5 Fr. 7 Fr.	40 cm 40 cm
05	Flexible Scissor Forcep	DM-S5-304-40 DM-S7-304-40	5 Fr. 7 Fr.	40 cm 40 cm
06	Semi Rigid Scissor Forcep	DM-R-S5-305-40 DM-R-S7-305-40	5 Fr. 7 Fr.	40 cm 40 cm



FORCEP FOR NEPHROSCOPY



S.NO.	DESCRIPTION	ART NO.	DIAMETER	WORKING LENGTH
01	Alligator Forcep for P.C.N.L.	DM-PC-25-306	2.5 mm	35-40 cm
		DM-PC-30-306	3 mm	35-40 cm
		DM-PC-35-306	3.5 mm	35-40 cm
02	Bi-prong Forcep for P.C.N.L.	DM-PC-25-307	2.5 mm	35-40 cm
		DM-PC-30-307	3 mm	35-40 cm
		DM-PC-35-307	3.5 mm	35-40 cm
03	Triprong Forcep for P.C.N.L.	DM-PC-25-308	2.5 mm	35-40 cm
		DM-PC-30-308	3 mm	35-40 cm
		DM-PC-35-308	3.5 mm	35-40 cm
04	Alligator Forcep for URS	DM-UR-40-309	4 Fr.	60 cm
		DM-UR-50-309	5 Fr.	60 cm
05	Triprong Forcep for URS	DM-UR-30-310	3 Fr.	60 cm
		DM-UR-40-310	4.2 Fr.	60 cm



ACCESSORIES

S.NO.	DESCRIPTION	ART NO.
01	Dumble Washer	DM-CYS-200-06
02	Patient Plate Cord	DM-CYS-200-08
03	Patient Plate long steel	DM-SL-200-09
04	Patient plate short steel	DM-SS-200-10
05	Pencil switch with cord & jack	DM-PS-200-11
06	TUR Adopter	DM-TU-200-12
07	Silicon Patient Plate with Cable	DM-SP-200-13

