


Laparoscopes Standard & HD



Benefits

- All endoscope are compatible to the fittings of all common light guides.
- Eyepieces are compatible to all common couplers.
- All HD Endoscopes are manufactured using Specialized High Definition Lens array.
- All metal parts are made of stainless steel.

Diameter	Length	Angle Of View	Type	Art. No
10 mm	310 mm	0°	Std.	PM-LS-10-0
10 mm	310 mm	0°	HD	PM-LS-10-0-HD
10 mm	310 mm	30°	Std.	PM-LS-10-30
10 mm	310 mm	30°	HD	PM-LS-10-30-HD
10 mm	300 mm	30°	HD	PM-LS-10-45
6.5 mm	310 mm	30°	Std.	PM-LS-65-30
5 mm	300 mm	0°	Std.	PM-LS-5-0
5 mm	300 mm	0°	HD	PM-LS-5-0-HD
5 mm	300 mm	30°	Std.	PM-LS-5-30
5 mm	300 mm	30°	HD	PM-LS-5-30-HD
3 mm	300 mm	0°	HD	PM-LS-3-0-HD
3 mm	300 mm	30°	HD	PM-LS-3-30-HD