



HU32 Reusable Uretero-Renoscope

Technical Specifications

Optical System	Pixels: 160K
	Field of view: 120°
	Depth of view: 3-100mm
	Illumination: optical fiber
Insertion Tube	Maximum outer diameter: 8.5fr
	Working length: 650mm
Instrument Channel	Working channel: 3.6fr
Bending Section	Angulation range: 285°
Function	Photo
	Video
	White Balance
	Storage: 16G
Connections	DVI
	S-Video
	CVBS

AN ENDOSCOPE SOLUTION COMPANY

**HU32 Reusable
Uretero-Renoscope
Always Available
Cost Effective**



Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan
Community, Yuanshan Street, Longgang District
Shenzhen, Guangdong Province, 518115, China
Tel: +86-755-2227 5833
Fax: +86-755-2227 5833
Email: sale@hugemed.net

www.hugemed.net

111-0022-00 Ver 1.0



HU32 Reusable Uretero-Renoscope

Hugemed all in one system,
always be ready for portable image
processor



285°
Excellent
Deflection

For easier navigation
of lower pole


- Challenging Anatomy
- Passive Bending



Flexible Sheath
3.6Fr working channel



160K High Resolution



Dual Buttons
For snap / Video capture



Self Locker
For reducing surgeon fatigue



True 1: 1
Torque ratio