



LED LIGHT SOURCE

LED-1500



ISO : 9001 : 2015
ISO : 13485 : 2016



Salient Features:-

- Focus coherent lens technique for uniform brightness
- Continuous variable light intensity 1-100 %
- LCD Display of light intensity and lamp life
- Low noise level enabling quiet operation
- Compact and Light weight construction
- Extremely high optical output of light
- Optional remote control facility for essential functions.
- Special heat filter
- Standby mode

Technical Specifications:-

- Type of LED : Large area LED chip, ROHS, REACH compliant for light
- Illumination : Equivalent to 300w xenon light
- Color Temperature : $\geq 6500K$
- Standard Connection : Storz type, other types on request
- Classification : Class 1 / Class A Medical device
- Ambient temperature : $+10^{\circ}C$ to $+40^{\circ}C$ during storage
- Dimensions (W x H x D) : 35 x 13 x 37cm
- Power Supply : 100 / 240 \pm 10% VAC, 50/60 Hz
- Fuses : 2 x 3 Amp, 250 VAC
- Atmospheric pressure : 86 kPa-106kPa

Manufactured by:-



Divine Meditech Pvt. Ltd.

Plot No.: 18, Phase III, Industrial Area,
Tahliwal Dist. Una, H.P. - 174301 (India)
Contact: +91(9818396505)
E-mail: info@divinemeditech.com
Website: www.divinemeditech.com

Member:-



Authorized Distributor

