





Cystoscope



Products Details

Modelno.		Diameter	Working Length
PM-CYS-4-0-HD	0°	4MM	302MM
PM-CYS-4-30-HD	30°	4MM	302MM

URETERORENOSCOPE 6/7.5 Fr.



Double layer water blocking design


- Stone therapy with laser, ultrasound, electrohydraulic or ballistic lithotripsy
- Rigid auxiliary instruments can be used through the straight instrument channel
- Atraumatic distal end with smooth round tip for patient safety and ease of insertion.

Specifications of URS 6/7.5 Fr.

- Distal tip of sheath : 6 Fr.
- Type : Autoclavable with angled eyepiece
- Working Chanel : 4 Fr.

ART NO.	Telescope	Working Length	Irrigation port
DM-URS-001	12°	430 mm	2 Nos.
DM-URS-005	12°	315 mm	2 Nos.
DM-URS-006	12°	150 mm	2 Nos.

URETERORENOSCOPE 4.5/6.5 Fr.




- Stone therapy with laser, ultrasound, electrohydraulic or ballistic lithotripsy
- Rigid auxiliary instruments can be used through the straight instrument channel
- Atraumatic distal end with smooth round tip for patient safety and ease of insertion.

Specifications of URS 4.5 Fr.

- Distal tip of sheath : 4.5 Fr.
- Type : Autoclavable with angled eyepiece
- Working Chanel : 3.3 Fr.

ART NO.	Telescope	Working Length	Irrigation port
DM-URS-002	5°	430 mm	2 Nos.
DM-URS-003	5°	315 mm	2 Nos.

URETERORENOSCOPE 8/9.8 Fr.




- Stone therapy with laser, ultrasound, electrohydraulic or ballistic lithotripsy
- Rigid auxiliary instruments can be used through the straight instrument channel
- Atraumatic distal end with smooth round tip for patient safety and ease of insertion.

Specifications of URS 8/9.8 Fr.

- Distal tip of sheath : 8 Fr.
- Type : Autoclavable with angled eyepiece
- Working Chanel : 5 Fr.

ART NO.	Telescope	Working Length	Irrigation port
DM-URS-004	5°	430 mm	2 Nos.
DM-URS-007	5°	315 mm	2 Nos.
DM-URS-008	5°	150 mm	2 Nos.

NEPHROSCOPE



Double layer water blocking design










- Atraumatic distal end with smooth round tip for patient safety and ease of insertion.
- Compatible with operative sheath with hollow obturator for continuous irrigation and suction.

Specifications of Nephroscope

- Type : Autoclavable with angled eyepiece

Product	Art No.	Telescope	Operating Sheath	Working Channel
Rigid Mini Nephroscope	DM-PCNL-002	12 Fr.	15 Fr.	6 Fr.
Rigid Nephroscope	DM-PCNL-001	18 Fr.	21 Fr.	10.5 Fr.
Rigid Nephroscope	DM-PCNL-003	21 Fr.	24 Fr.	10.5 Fr.

TURP SET

	Description	Item No.
	Telescope 4mm 0° 4mm 30°	PM-CYS-4-0-HD PM-CYS-4-30-HD
	Working Element Monopolar Working Element Bipolar	PM-URA02799 PM-URA02799-SL
	Resectoscopy Sheath Outer Sheath 26 Fr. Inner Sheath 24 Fr.	PM-URD02798 PM-URD02797
	Resectoscopy Obturator Standard Obturator Visual Obturator	PM-URO02499 PM-URO02497
	Cysto- Urethroscope sheath 20 Fr.	PM-UCS02099
	Cysto- Urethroscope Obturator	PM-UCO02099
	Cystoscopy Bridge Single Channel	PM-UCU00198
	Urethrotome Sheath 21 Fr.	PM-UUS02199
	Urethrotome Obturator	PM-UUO02199